Lesson Element

Explaining Observations: rates of reaction of

halogenoalkanes

Instructions and answers for teachers

These instructions should accompany the OCR resource 'Explaining Observations: rates of reaction of halogenoalkanes' activity which supports OCR A Level Chemistry OCR A Level Chemistry B (Salters).



The Activity:

The simple, small-scale practical activity described here provides students with the opportunity to explore the precipitation reactions between solutions of halide salts and silver nitrate before going on to investigate the rates of hydrolysis of 1-chlorobutane, 1-bromobutane and 1-iodobutane.



This activity offers an opportunity for English skills development.



This activity offers an opportunity for maths skills development.

Associated materials:

'Explaining Observations: rates of reaction of halogenoalkanes' activity.





Overview

Understanding the properties and reactions of organic compounds is important if chemists are to develop safe and effective new materials using environmentally-friendly procedures.

The halogenoalkanes (also known as haloalkanes or alkyl halides) are found in many important industrial products as well as being useful (ie reactive) intermediates for conversion into other substances.

The effect of different halogen atoms (X) on the reactivity of some halogenoalkanes, R-X, can be investigated by comparing the rates of hydrolysis of 1-chlorobutane, 1-bromobutane and 1-iodobutane

$$.R-X + H_2O \rightarrow R-OH + H^+ + X^-$$

The reaction is carried out in in the presence of silver ions, Ag^+ . These react with the halide ions produced, X⁻, to form a precipitate, AgX.

$$Ag^{+}(aq) + X^{-}(aq) \rightarrow AgX(s)$$

The time taken to form the precipitate gives an indication of how quickly the halogen atom, X, in the halogenoalkane is replaced by an –OH group.

Students' common misconceptions

In order to suggest, develop and evaluate explanations for observed differences in the rates of reaction students need to draw on their understanding of (i) ionic and covalent compounds, (ii) electronegativity and bond polarity and (iii) bond enthalpy (bond strength).

However, even though they may have previously used silver nitrate solution to test for the presence of halide ions, some students do not fully appreciate the significance of the silver halide precipitate and the clues its formation provides about the hydrolysis reactions of halogenoalkanes.

There is evidence that some post-16 students are unable to explain what is happening at the level of particles when precipitation occurs and that this may be related to their lack of understanding of what happens during dissolving and their ideas about the nature of ionic bonding.

[Adapted from: Chemical misconceptions – prevention, diagnosis and cure. Volume II by Keith Taber. Royal Society of Chemistry 2002]

Teachers may like to use the diagnostic probe activities (see links below) with the practical activity described here to explore their students' existing understanding of some of the key ideas involved. Teacher notes and expected answers are included with the probes. Some teachers may also be





interested in reading the more in-depth background research findings relating to the use of these probes with groups of pre-16 and post-16 students.

Diagnostic probes (classroom resources)

- Precipitation: http://www.rsc.org/learn-chemistry/resource/res00001096/precipitation
- Ionic bonding: http://www.rsc.org/learn-chemistry/resource/res00001095/ionic-bonding

Theoretical background (research findings)

- Students' ideas about precipitation processes: http://www.rsc.org/learnchemistry/resource/res00001141/chemical-reactions
- Students' ideas about ionic bonding: http://www.rsc.org/learnchemistry/resource/res00001140/chemical-bonding

Teacher and Technician notes

Practical requirements

Each student (or pair of students) will require:

- Eye protection
- Beaker, 100 cm³
- Small test-tubes or combustion tubes, 75 mm x 10 mm (4)
- Student worksheet plus a laminated drop-chemistry/data card (or put inside a plastic pocket)
- Plastic dropper pipettes, 1 cm³ (2)
- Thermometer, 0-110 °C

Students will also require access to:

- Hot water (use a kettle to avoid sources of ignition/naked flames)
- Silver nitrate solution, 0.05 mol dm⁻³ [Solution will stain skin]
- Ethanol or IDA/IMS [Highly flammable, Harmful if swallowed]
- Dropper bottles (or equivalent) containing:
- Sodium (or potassium) chloride, bromide and iodide solutions, each at 0.1 mol dm⁻³
- 1-chlorobutane [Highly flammable]
- 1-bromobutane [Highly flammable, Irritating to eyes, respiratory system and skin]
- 1-iodobutane [Flammable, Harmful by inhalation and in contact with skin]





Task 1 Investigating the precipitation reaction with halide salts

Wear eye protection.

Carry out the following reactions on the laminated card circles. Note the colours of the precipitates formed.

- Circle A: 1 drop of chloride solution + 2 drops of silver nitrate solution.
- Circle **B**: 1 drop of bromide solution + 2 drops of silver nitrate solution
- Circle C: 1 drop of iodide solution + 2 drops of silver nitrate solution

http://www.rsc.org/learn-chemistry/resource/res00001096/precipitation

The notes accompanying the Precipitation diagnostic probe activities (see hyperlink) provide some possible answers that post-16 students could be expected to suggest (see activities A reaction to form silver chloride and A precipitation reaction).

Students may offer different explanations depending on which topics they have already studied (eg, explanations in terms of energy changes in solution).

Students could just be provided with the diagrams from the probes to use in their explanations.

Task 2 Investigating the rate of hydrolysis of halogenoalkenes

Students should recognise that this time a precipitate does not immediately occur because the halogen atoms are covalently bonded in the halogenobutane and are not already present in solution as solvated halide ions. As the hydrolysis (a nucleophilic substitution reaction) progresses, halide ions are gradually produced.

In trials at about 40 °C the iodobutane produced a precipitate within a few seconds - this gradually became thicker. The bromobutane took around 1-2 minutes to produce the first signs of a precipitate and the chlorobutane produced a faint precipitate after about 5-7 minutes.





The order in which precipitates appear is explained by the trend in bond enthalpies rather than bond polarities (ie bond has lowest bond enthalpy, bond breaks most quickly so produces precipitate first).

Drawing the ideas together, students should be able to offer some logical ideas for a possible mechanism for the hydrolysis reaction.

Teachers may like to provide students with an opportunity to model the nucleophilic substitution mechanism. A set of cards and additional information can be found in Activity 55 on the Royal Society of Chemistry's *Assessment for learning chemistry* website,

<u>http://www.rsc.org/education/teachers/resources/aflchem/</u>. Some students have produced effective animations of the mechanism using various software tools.

We'd like to know your view on the resources we produce. By clicking on '<u>Like</u>' or '<u>Dislike</u>' you can help us to ensure that our resources work for you. When the email template pops up please add additional comments if you wish and then just click 'Send'. Thank you. If you do not currently offer this OCR qualification but would like to do so, please complete the Expression of Interest Form which can be found here: <u>www.ocr.org.uk/expression-of-interest</u>

OCR Resources: the small print

OCR's resources are provided to support the teaching of OCR specifications, but in no way constitute an endorsed teaching method that is required by the Board, and the decision to use them lies with the individual teacher. Whilst every effort is made to ensure the accuracy of the content, OCR cannot be held responsible for any errors or omissions within these resources. We update our resources on a regular basis, so please check the OCR website to ensure you have the most up to date version.

© OCR 2015 - This resource may be freely copied and distributed, as long as the OCR logo and this message remain intact and OCR is acknowledged as the originator of this work.

OCR acknowledges the use of the following content: Maths and English icons: AirOne/Shutterstock.com Please get in touch if you want to discuss the accessibility of resources we offer to support delivery of our gualifications: resources.feedback@ocr.org.uk





